



MONet XCT Working Group Overview of XCT workflow and data


Tamas Varga

Billy Petersen

tamas.varga@pnnl.gov

billy.petersen@pnnl.gov

Sep 19, 2024




- Brief intro to X-ray Computed Tomography (XCT)
- MONet XCT workflow
- MONet XCT data types
- Other data output formats/options, software used

What is XCT?


- Non-destructive x-ray imaging to examine/visualize the internal microstructure of objects
- Output: a series of radiographs computationally reconstructed as 3D volume (CT)
- Materials are differentiated in image data based on their electron density (x-ray absorption); difference in x-ray absorption results in different gray level (color) intensity in the image



How does our x-ray CT work?



Cone Beam Configuration




http://serc.carleton.edu/images/research_education/geochemsheets/techniques/_1172960774.png


<https://www.emsl.pnnl.gov/science/instruments-resources/x-ray-computed-tomography>

A few key terms


- Voxels (size), volume data



- Slices/stacks




- Segmentation




- Gray-scale vs. RGB




- Mesh



- Radiograph




MONet core processing workflow




Scanning protocol


Low resolution XCT
scans collected on
whole core in 2
tiles
(82.5 microns voxel
size)



High resolution XCT
scans collected on
top and bottom
sections
(38.2 microns voxel
size)



Whole core (merged low-resolution data from top and bottom sections)




Tomographic image of
whole core in XZ
plane view


The two lower-resolution (82.5-micron voxel size, estimated spatial resolution 165 microns) scans are used to get the merged image of the whole 1-ft tall core. This image has been stitched together from two tiles as suggested in the title.

Scan conditions:

X-ray power: 105 kV / 325 μ A
Projections: 1600
Frames/projection: 2
Exposure time: 500 ms
Filter: 0.5-mm Cu foil
Voxel size: 82.5 μ m



Higher-resolution scans




Tomographic images
of top and bottom
sections of core in XZ
plane view

The two higher-resolution (38.2-micron voxel size, estimated spatial resolution 76.4 microns) scans are used to perform full pore network characterization. These images cannot be stitched together into one.

Scan conditions:



X-ray power: 105 kV / 325 μ A
Projections: 3142
Frames/projection: 2
Exposure time: 708ms
Filter: 0.5-mm Cu foil
Voxel size: 38.2 μ m




- 📁 60880 Libohova_MONet Core 4 A top zoom-files
- 📄 60880 Libohova_MONet Core 4 A top zoom.hx
- 🗄️ 60880 Libohova_MONet Core 4 A top zoom_Data
- 60880 Libohova_MONet Core 4 A top zoom_Pores connected
- 60880 Libohova_MONet Core 4 A top zoom_Pores PNM
- 60880 Libohova_MONet Core 4 A top zoom_Pores Unconnected
- 60880 Libohova_MONet Core 4 A top zoom_Pores
- 60880 Libohova_MONet Core 4 A top zoom_Raw

3D image reconstruction


- Create 3D volume data from the individual projections
- We use CT Pro 3D (X-Tek software suite by Nikon)




Higher-resolution scan (with 38.2-micron voxel size) on the top section of core




Reconstructed, cropped volume




Visualizing all pores



Connected pores only



Isolated pores only




- Total Porosity: 10.97 %
- Connected Pores Percentage: 89.75 %
- Wet Bulk Density: $\approx 2.05 \frac{g}{cm^3}$


- Mean Pore size: 0.06448 mm
- Min Pore size: 0.047327 mm
- Max Pore size: 29.63071 mm
- Median Pore size: 0.047338 mm
- Variance Pore size: 0.001942 mm
- Mean Pore Volume: 0.005736 mm³

Cropped volume size analyzed here:
50x45x60 mm


Pore network modeling (PNM) results on higher-resolution data from top section of core




Reconstructed, cropped volume



Visualizing all pores



Connected pores



PNM

Cropped volume:

50x45x60 mm

Input pressure: 103098 Pa,
Output pressure: 102833 Pa,
Fluid viscosity: 0.001 Pa.s

Pore Network Model in Z-direction

- Absolute permeability: $109.69 \mu\text{m}^2$
- Total flow rate: 1102.00 mm^3
- Tortuosity: 1.80

Pore Network Model in Y-direction

- Absolute permeability: $2882.55 \mu\text{m}^2$
- Total flow rate: 50706.88 mm^3
- Tortuosity: 1.52







Pore Network Model in X-direction

- Absolute permeability: $379.04 \mu\text{m}^2$
- Total flow rate: 5450.30 mm^3
- Tortuosity: 1.76



Data types for MONet

- Image data:




New: binary stacks!

- Numerical data:

- Pore size
- Pore volume fraction
- Pore connectivity
- Pore size distribution
- Bulk density
- Permeability
- Flow rate
- Tortuosity

Our methodology to estimate bulk density from CT data



Estimate of bulk density is based on formula from Tanaka *et al.*, *Earth Planets Space*, **63**, 103–110, 2011:

$$HU = \frac{1000(\mu - \mu_w)}{\mu_w},$$

where HU means Hounsfield Units, μ is the grey value of a given voxel and μ_w is the grey value of water for our scan. μ_w was obtained from finding the average grey value of water in a previous scan with similar conditions. Then we used *HUs* to estimate the bulk density for a given pixel through the formula

$$\rho = \frac{0.677}{1000} HU + 1.0756.$$

Finally, we averaged the bulk densities for all pixels that are soil or pore (excluding wall porosity) to find the overall average bulk density of the scan.

Wet bulk density: $\approx 1.50 \frac{g}{cm^3}$


Data processing

- Importing reconstructed data into Avizo (Thermo Fisher)
- Image cropping
- Filtering
- Segmentation
- Pore volume/size analysis
- Pore network analysis
- Numerical output
- Visualization

Other data output options

- **Image data:**

- Animations
- Mesh
- STL file (3D printing)
- ...



- **Numerical data:**

- Pore throats
- Pore coordination numbers

Viewing raw and segmented data (ImageJ/Fiji)

Google search results for "download fiji".

Search bar: download fiji

Filter buttons: Free, For windows 7, Ringtone, Videos, Song, ImageJ, Images, Mac, Radio

Results:

- ImageJ Wiki**
https://imagej.net › software › fiji › downloads
- Fiji Downloads**
Fiji is a distribution of ImageJ which includes many useful plugins contributed by the community.
~ [Download Fiji](#) for your OS ~. Windows 64-bit, imagej.net (USA) ...
[Archive · Life-Line Fiji versions](#)
- Fiji**
https://fiji.sc
- Fiji: ImageJ, with "Batteries Included"**
Fiji is an image processing package — a "batteries-included" distribution of ImageJ, bundling many plugins which facilitate scientific image analysis. [Download ...](#)
- ImageJ Wiki**
https://imagej.net › software › fiji
- Fiji**
Downloads. ~ [Download Fiji](#) for your OS ~. Windows 64-bit, imagej.net (USA), micron.ox.ac.uk (European mirror). Windows 32-bit, imagej.net (USA), micron ...

A blue arrow points to the "Fiji: ImageJ, with 'Batteries Included'" result.



Install Fiji

Fiji is an image processing package — a "batteries-included" distribution of [ImageJ](#), bundling many plugins which facilitate scientific image analysis.

[Download for Windows \(64-bit\)](#)

[More Downloads](#) [Cite](#) [Contribute](#)

Why Fiji?



Easy to Use





Powerful




Free & Open Source

Install Fiji – cont.


- Copy Fiji.app folder from zip file to C drive
- Run application file “ImageJ-win64”
- Allow for Java installation as offered (if you do not already run Java 1.8)
- Once app started, run updater
- You can pin this application file to taskbar for easy access





Importing raw scan data (.vol file)




Importing data – cont.





Enhancing image contrast






Enhancing image contrast – cont.





Filtering



Importing segmented, binary image stack





Importing stack – cont.




Creating STL file from segmented, binary image data


- Step-by-step process from binary stack (previous slide) to .stl in ImageJ



Creating STL file from segmented, binary image data – cont.




Creating STL files from segmented, binary image data – cont.



3A_bot_all_pores_binary_stack	5/20/2024 9:09 AM	File folder	
3A_bot_conn_pores_binary_stack	5/20/2024 9:09 AM	File folder	
3A_bot_unconn_pores_binary_stack	5/18/2024 11:36 AM	File folder	
61111 Scott MONet Core 3 A bot zoom-files	5/6/2024 3:55 PM	File folder	
61111 Scott MONet Core 3 A bot zoom PNM.tif	5/6/2024 3:57 PM	TIFF image	233 KB
61111 Scott MONet Core 3 A bot zoom PORES.1241x1201x1576 8us.raw	5/6/2024 4:03 PM	RAW File	2,293,882 KB
61111 Scott MONet Core 3 A bot zoom PORES.CONNECTED.tif	5/6/2024 4:01 PM	TIFF image	429 KB
61111 Scott MONet Core 3 A bot zoom PORES.unCONNECTED.tif	5/6/2024 4:02 PM	TIFF image	388 KB
61111 Scott MONet Core 3 A bot zoom PORES.Axis-Connectivity.1241x1201x1576 8us.raw	5/6/2024 4:08 PM	RAW File	2,293,882 KB
61111 Scott MONet Core 3 A bot zoom PORES.tif	5/6/2024 3:59 PM	TIFF image	466 KB
61111 Scott MONet Core 3 A bot zoom PORES.xor.1241x1201x1576 8us.raw	5/6/2024 4:09 PM	RAW File	2,293,882 KB
61111 Scott MONet Core 3 A bot zoom RAW.tif	5/6/2024 3:58 PM	TIFF image	427 KB
61111 Scott MONet Core 3 A bot zoom.hx	5/6/2024 3:56 PM	Avizo Project	55 KB
61111 Scott MONet Core 3 A bot zoom.xls	5/6/2024 4:27 PM	Microsoft Excel 97...	6,106 KB
Conn_porees.stl	9/13/2024 11:50 AM	Certificate Trust List	43,479 KB


Viewing STL files

- STL Viewer: <https://www.viewstl.com/>



Viewing STL files

- STL Viewer: <https://www.viewstl.com/>



- Avizo: <https://www.thermofisher.com/us/en/home/electron-microscopy/products/software-em-3d-vis/avizo-software.html>
- ImageJ/Fiji: <https://imagej.nih.gov/ij/download.html> and <https://fiji.sc/>
- ParaView: <https://www.paraview.org/download/>
- Modeling:
 - OpenPNM - Pore Network Modeling Framework in Python:
 - <https://openpnm.org/>
 - <https://github.com/PMEAL/OpenPNM>
 - PoresPy - Quantitative Image Analysis of Porous Materials:
 - <https://porespy.org/>
 - PFLOTRAN – A massively parallel subsurface flow and reactive transport code:
 - <https://www.pfotran.org/>
 - OpenFOAM - A free, open source CFD software:
 - <https://www.openfoam.com/>

- Visual Graphics Studio, Dragonfly (license)
- ImageJ/Fiji, ParaView (free)
- Modeling:
 - OpenPNM:
 - <https://openpnm.org/>
 - <https://github.com/PMEAL/OpenPNM>
 - PoresPy:
 - <https://porespy.org/>
 - PFLOTRAN:
 - <https://www.pfotran.org/>
 - OpenFOAM:
 - <https://www.openfoam.com/>



Thank you!



MONet

[emsl.pnnl.gov
/MONet](http://emsl.pnnl.gov/MONet)

- DOE-BER
- MONet Team



J. Bargar



E. Graham



S. Leichty



O. Qafoku



A. Townsend



Y. Corilo



S. Karra



Y. Song



M. Rockhold



T. Wietsma



J. Toyoda



B. Petersen



M. Mudunuru



E. Elof-Fadrosch



K. Thibault

... and many more

Questions?

